Viewer Lite

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User Manual

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Viewer Lite

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8300234036 Rev. A

Esaote S.p.A.
Genova

Edizione

Settembre 2008

Garanzia

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September 2008

Guarantee

The information contained in this documentation is the property of Esaote S.p.A. All rights reserved. Reproduction and distribution in any form is strictly prohibited. All rights reserved

Translations of this document are strictly prohibited without prior authorization from Esaote S.p.A..

The information contained in this documentation is subject to change without prior notice.

Trade Marks

DICOM3® is a registered trademark of National Electrical Manufacturers Association

Windows® is a registered trademark of Microsoft Corporation

All names are property of the respective owners and are used exclusively for identification purposes.

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Main authorized distributors

HOLOGIC INC
35 Crosby Drive, Bedford
MA 01730, U.S.A.
Tel. +1 781 9997300

ITISA BIOMEDICA
Avenida Camino de lo Cortao 34, nave 7
28700 - San Sebastian de los Reyes, Spain
Tel. +34 91 6592202 Fax +34 91 6512236

SCANEX Medical System AB
La Cours Gata 4
S-252 31 Helsinborg, Sweden
Tel. +46 42 373400 Fax +46 42 242472

SCANEX Medical System A/S
Solheimveien 112
N 1473 Skarer, Norway
Tel. +47 67 921150 Fax +47 67 921151

OY SCANEX Medical System AB
Solheimveien 9
P.O. Box 55
FIN 02361 ESPOO, Finland
Tel. +358 9 8886007 Fax +358 9 88170512

SCANEX Medical System A/S
Forskningscentret Venlighedsvej 6
DK - 2970 Horsholm, Denmark
Tel. +45 45 766661 Fax +45 45 766697

MEDISON ECONET CO. LTD.
Daeyoung DB 4F, 9-1 Samsung-dong, Gangnam-Gu,
135-092 Seoul, Korea
Tel. +82 -2 -34161904 Fax +82 -2 -34161969

IBERDATA EQUIPAMENTO sa.
R. Dr. Ricardo Jorge, 4
Venda Nova 2700 Amadora, Portugal
Tel. +351 21 4768176 Fax +351 21 4768179

ORTOPEDIA HISA
Gutenberg No. 138 Col. Anzures, Mexico D.F.
Tel. +52 55319495 Fax +52 55319989

BRACCO UK LTD
Bracco House HP10 OHH
Mercury PK, Wycombe LN, Wooburn Green
High Wycombe Bucks HP10 - GB
Tel. +44 1628851515 Fax +44 1628851505

BIOLAB (Aust) Pty Ltd
2 Clayton Road Clayton VIC 3168 Locked Bag 24
Mulgrave VIC 3170- Australia
Tel. +61 392634300 Fax +61 395487060

ATTIEH MEDICO
P.O. Box 116105
Jeddah 21391 - Saudi Arabia
Tel. +966 26972550 Fax +966 26821582

TAIYO DENSHI CO. LTD.
140-2 Nishibaramae Oya Ogawara Machi Sibata-Gun
Miyagi Ken - Japan
Tel. +81 224521846 Fax +81 224522715

SUPERVISION SYSTEM

The Viewer Lite equipment is subject to a supervision system (post-sales supervision), which ESAOTE S.p.A., all associates and authorized distributors apply to products issued onto the market, in relation to real or potential hazards that may arise for the patient or operator during normal use of the equipment, to ensure optimal solutions in the most efficient and prompt manner possible.

Therefore, in the event of malfunctions, defective performance of the equipment, or inadequate instructions, which may constitute a hazard to the patient or operator, the user must notify ESAOTE S.p.A., associate company or authorized distributor in writing, providing the information as specified in the form below. The data of Viewer Lite can be viewed in the window "About Viewer Lite".

On receipt of the notification Esaote S.p.A. will immediately activate the process of examination and resolve the non-conformity that has been reported.

ACCIDENT REPORT FORM

TO: ESAOTE S.p.A.
Quality Assurance Service
Via Siffredi 58
16153 Genoa, Italy
[or associate company]
[or authorized distributor]

Device Name

ESAOTE

Code (REF)

Serial Number (S/N).....

Description of accident or potential accident

.....

Comments or suggestions

.....

Person / department to contact

Address

Telephone Fax

Date

Signature

USER LICENSE CONDITIONS FOR VIEWER LITE SOFTWARE

Warning

Carefully read the terms and conditions below before using the software.

Use of the software implies acceptance of the terms and conditions stated below.

PROPERTY

Viewer Lite includes software property of Esaote S.p.A. and/or software granted by license to Esaote S.p.A. by one or more suppliers ("Software suppliers"). These software products ("SOFTWARE"), including all relative media, printed material and "online" documentation or electronic material, are protected by law and international agreements on Intellectual Property. The SOFTWARE is not sold, but granted on license. The SOFTWARE, and, similarly, any copyrights and all industrial and intellectual ownership rights are and shall remain the exclusive propriety of Esaote S.p.A. or its third party suppliers.

The user does not acquire rights to the SOFTWARE, with the exception of the granted user license.

USER LICENSE AND LIMITATIONS

With this license, Esaote S.p.A. grants the final user the right to use the software on one computer only.

In any event the user may not make unauthorized copies and/or reproductions of the SOFTWARE or parts of the latter, including the documentation enclosed.

On the basis of the above, and given that the SOFTWARE is not protected against copying, one copy only may be made for safety purposes (backup).

The user may not hire, lease or loan out the SOFTWARE, but may transfer rights permanently according to the conditions herein, provided that all copies of the SOFTWARE and written material are transferred and the new owner accepts all conditions. Transfer of the product must always include the latest update and all previous versions.

The user may not convert, decode, disassemble or modify the SOFTWARE.

The user may not remove, tamper with or alter the copyrights, trademarks of other proprietary rights affixed or enclosed in the SOFTWARE.

The user may not publish data or information that compares performance of the said SOFTWARE with that written by others.

INDEX

Chapter 1 - Essential Prescribing Information

Intended Use	1
Manual Requirements	2
Safety Information: Conventions	2
Safe Use - Responsibility.	2
PC Requirements	3
Compatibility with some Internet browsers for Windows Vista . . .	3
Equivalence between European and international standards . . .	4

Chapter 2 - Viewing and saving images

Workspace	1
Displaying images	1
Selecting a series of images	1
Displaying the required series	2
Selecting the image display format	3
Displaying more than one series at a time.	6
Arranging windows in the image viewer	6
Image identification information	8
Image layout settings.	10
Browsing between images.	11
Displaying previous and next images	11
Displaying previous and next series	12
Image properties	13
Saving images in non-diagnostic formats	13

Chapter 3 - Image processing functions

Enlarging images	1
Changing image brightness and contrast	2

Image pan	3
Changing image orientation	4
Image Clips.	5
Statistics.	7
Measurements.	9
Measuring the distance between pairs of points	9
Measuring the angle between two segments	10
Inserting text and direction arrows.	11
Deleting measurements taken	12
Deleting all changes made to an image	13

Chapter 4 - User Interface management

Showing / hiding tool bars	1
Window management.	1
Closing a window	3
About VLite	4

CHAPTER 1

Essential Prescribing Information

.....

The software Viewer Lite is produced by Esaote S.p.A.
Via Siffredi 58, 16153 Genoa, Italy
Telephone +39 10 65471 FAX +39 10 6547275

The software package Viewer Lite is an option included in the software version installed on Esaote magnetic resonance equipment.

The software version is specified in the window "about VLite" (see chapter 5 page 4).

Intended Use

Viewer Lite is a software used to view and save digital images acquired by Esaote magnetic resonance devices.

Warning

Viewer Lite does not replace the skills and judgement of a qualified doctor who has been suitably and adequately trained in the software functionalities, and must be used in strict observance of the contents of this manual.

Warning

**The user must be aware of the limitations regarding the accuracy and correctness of data viewed and/or saved using Viewer Lite.
The quality of the data depends on the correctness of input data, interaction of the user with the data, the quality, characteristics and settings of the PC and the monitor on which the information is viewed.**

Warning

Users are not permitted to use the software in breach of the limits as established by current local and/or federal and/or state legislation.

Manual Requirements

This User Manual, code 8300234036 Rev. A, describes use of the Viewer Lite software and provides all instructions required to ensure safe and efficient use of the product.

The original language of all the E Viewer Lite Software User Manuals is Italian. Esaote S.p.A. has translated these manuals into English and German.

Safety Information: Conventions

The following labelling conventions are used in this manual to identify safety related information.

WARNING

A warning message alerts the reader to a situation which, if not avoided, could result in death or serious injury to the patient, the user or others. It may also describe potentially serious adverse reactions and safety hazards.

NOTE

A note informs the reader of information concerning the installation, operation or maintenance which is important, but not hazard-related.

Safe Use - Responsibility

Esaote S.p.A. assumes no responsibility for personal injury or property damage arising from the User's failure to observe the instructions in the User Manual (and any attachments), which must be available - to operation personnel - intact and in fully legible conditions.

Users must carefully read this User Manual, containing a detailed description of use of the various software functions.

If the user is uncertain with regard to the correct and safe use of Viewer Lite, in the event of malfunctions or deterioration of the software itself, or when information in this manual is incorrect or incomplete, proceed as follows:

- avoid use of the software
- contact Esaote S.p.A. authorized service personnel for assistance.

PC Requirements

Viewer Lite can be installed on any commercial PC that meets the following requirements:

- Operating System Windows 2000 or Windows XP Professional
- Microsoft Service Pack 4 for Windows 2000
- Internet Explorer 6.0 SP1 (minimum)
- Adobe Acrobat Reader 5.0 or later version
- Graphics card with 16MB VRAM
- 256 MB RAM
- CD/DVD drive and/or optical-magnetic disk drive.
- Color monitor, minimum requirements: 17", resolution 1280 × 1024, 24 bit, VGA = 60 Hz

Warning

The PC and monitor must comply with the safety standard EN 60601-1:1990 and subsequent amendments, or EN 60950-1:2001 and subsequent amendments. If the devices comply with the safety standard EN 60950-1 they should be installed away from the patient's area, as foreseen by the safety standard EN 60601-1-1:2001 and subsequent amendments. The PC and monitor should also comply with safety standard 60601-1-2:2001 and subsequent amendments or the electromagnetic compatibility provisions of 89/336/EEC Directive and subsequent amendments.

Warning

If the PC configuration is not equal to or greater than the minimum recommended requirements, correct operation of the software cannot be guaranteed.

Compatibility with some Internet browsers for Windows Vista

Using a program for browsing Internet is necessary only if the images in JPEG format were exported onto CD/DVD.

The link to the Viewer lite on the web page "index.htm" can not be used due to possible incompatibilities between the internet browser used and the Viewer Lite.

Equivalence between European and international standards

European Standards	International Standards
EN 60601-1:1990 Electromedical equipment - Part 1: General safety standards	IEC 60601-1:1988 Electromedical equipment - Part 1: General safety standards
EN 60601-1-1:2001 Electromedical equipment - Part 1: General safety standards – Collateral standard: Safety requirements for electromedical systems	IEC 60601-1-1:2000 Electromedical equipment - Part 1: General safety standards – Collateral standard: Safety requirements for electromedical systems
EN 60601-1-2:2001 Electromedical equipment - Part 1: General safety standards – Collateral standard: Electromagnetic compatibility – Requirements and tests	IEC 60601-1-2:2001 Electromedical equipment - Part 1: General safety standards – Collateral standard: Electromagnetic compatibility – Requirements and tests
EN 60950-1:2001 Information technology equipment – Safety – Part 1: General Requirements	IEC 60950-1:2001 Information technology equipment – Safety – Part 1: General Requirements

CHAPTER 2

Viewing and saving images

.....

The Viewer Lite lay-out is made up of the work window containing the patients and studies - on the left - and the display window, which loads and shows images related to series selected from the workspace.

Workspace

To open/close the window containing the workspace, from the main menu select **View** → **Workspace**

Displaying images

The viewer can only display series of images stored on the removable media, whether optical disk or CD/DVD.

Selecting a series of images

To select a series of images in this way, display images in this way, proceed as follows:

- 1 Insert the optical disk or CD/DVD in the relative drive
- 2 If the workspace on the left of the image viewer is not open, from the main menu select **View** → **Workspace**
- 3 The list is sorted in descending order, first by the patients in the directory DICOM® (DICOMDIR), followed by the list of any images that have not been saved on the media according to the DICOM3®, standard.

The DICOMDIR is organized in a tree structure of patients/studies/series.

Each patient folder contains the last name and name of the latter and clicking on + next the patient icon displays a list in cascade of the relative studies.

Each patient folder contains the last name and name of the latter and clicking on + next the patient icon displays a list in cascade of the relative studies.

Each series folder contains the time of execution of the series and information on the type of sequence acquired; including slice orientation.

- 4 select the required series and double click with the left mouse button on the line containing the specific series to open in the image viewer.


fig. 2.1 - Workspace

Note that the image viewer can open magnetic resonance images only. If an attempt is made to open different types of images (for example ultra-sound), a dialog box tells the user that this type of image cannot be opened.

Displaying the required series

Once the series to be displayed has been selected and opened according to the procedure described in the above paragraph, the first image of the series is displayed in a window to the right of the workspace.

The image series viewer window is divided into three areas:

- the first, large and located at the top, contains the image/s loaded as described above.

- the second, located below or to the right of the above, contains the three Scout images related to the loaded series (except when the loaded series is the Scout itself, in which case nothing is displayed).

The lines displayed on the three Scout images correspond to the intersection of each, with the levels making up the package used to acquire the sequence.

These are red, except for the one corresponding to the level that cuts the image in the currently displayed series, which is yellow.

The user can show/hide the area containing the three Scout images related to a loaded series, by following one of the procedures below:

- 1 METHOD 1, from the main menu, select **View** → **Scout Area**
- 2 METHOD 2, click on the corresponding icon on the tool bar displayed alongside.

The selection is not only linked to the current series but remains active until a new modification is made by the user.

Selecting the image display format

Images of the series may be displayed in three ways: Single image on full screen, 2×2 images and "custom" ($n \times m$) images). The user can change from one configuration to another as described below:

to view images in format 1

- 1 METHOD 1, from the main menu, select **View** → **1 × 1**
- 2 METHOD 2, click on the corresponding icon on the tool bar alongside


fig. 2.2 - Displaying a single image

to display images in 2×2 format

- 1 METHOD 1, from the main menu, select **View** → **2 × 2**
- 2 METHOD 2, click on the corresponding icon on the tool bar alongside


fig. 2.3 - 2 × 2 display

☞ to display images in custom format

- 1 METHOD 1, click on the corresponding icon on the tool bar displayed alongside. This displays a submenu to enable, by means of the mouse, the user to customize the display by choosing the number of rows and columns. Simply complete part of the table shown, which will become blue
- 2 METHOD 2, select **View** → **Custom...** from the main menu to automatically display a window like the one shown in the following figure, for the selection of the pre-defined format, by choosing the number of rows and columns.

fig. 2.4 - "Custom" display mode


fig. 2.5 - Custom format setup window

To change the number of rows and columns, click with the left mouse button on the ↑ and ↓ arrows located to the side of the boxes **Columns** and **Rows**, which enable the user respectively to increase or decrease the numbers of columns and the number of rows.

Alternatively, click with the left mouse button inside the field containing the value to be modified and enter the required number.

Displaying more than one series at a time

It is also possible to load different series belonging to the same patient study, different studies from the same patient or also different patients: the selected series are loaded in the image viewer as described below.

When the first series is loaded, as described above, a full screen window is displayed; when the second series is opened, a new window with the series is displayed on screen according to the window layout settings (see sub-paragraph below).

Arranging windows in the image viewer

When more than one series is opened at a time, the system, by default, places the relative windows side by side inside the image viewer: the second series opened is placed horizontally next to the first, which is redimensioned, after which if other series are opened, the display follows the tile display as adopted in the Windows®, operating system, and shown in the following figure.

fig. 2.6 - "Tile" arrangement of windows

If the user changes the arrangement of a window or its size, the viewer no longer automatically manages window positioning. The windows can be re-arranged by default at any time:

select **Window** → **Auto Arrange**, from the main menu or click on the relevant icon on the toolbar.

The windows containing new series to be displayed can also be arranged in cascade:

select **Window** → **Cascade**, from the main menu or click on the relevant icon on the toolbar displayed alongside.


fig. 2.7 - Displaying different series in cascade

To return to the automatic tile arrangement of windows:

select **Window** → **Auto Arrange**, from the main menu or click on the relevant icon on the toolbar.

With both the tile and cascade display, any window can be maximized to full screen as required.

To maximize a window, click with the left mouse button on the window's maximize/minimize button, or double click on the title bar.

To minimize a window, click with the left mouse button on the window's maximize/minimize button, or select the command **Window** → **Arrange all**.

Image identification information

The display information required to identify the image differs according to the display format.

When the display format is reduced, the system first reduces the font of all displayed information automatically, and then decreases the amount of information to improve image clarity. This is to accelerate and simplify interpretation of the image.

If an image is displayed full screen, the information will be shown as illustrated in the figure below.

Sequence name		Dicom		Last and first name		
TR / TE / NEX		ID		Birth date		Sex
FOV		District		Laterality		
Matrix						
Scan Time						
FA		Examination date				
TI		Examination hour				
# Echo						
# Slice						
Thickness						
Distance		Dicom		Dicom		
Hamming Filter						
Quality Factor						
Center						
Width						
Zoom						
System's name		Dicom			Institute	
Image					Department	

fig. 2.8 - Displaying information in "full screen" format

The information column to the left of the image displays:

- type of sequence
- repetition time in ms - echo time in ms - number of excitations
- acquisition matrix
- field of view (mm × mm)
- scan duration in minutes:seconds
- flip angle in degrees
- interpulse time in ms
- echo / total number of echoes
- slice / total number of slices
- thickness in mm
- distance between current slice and magnet isocenter in cm
- contrast bolus (when used)
- Hamming filter level used
- Quality factor expressed as a percentage
- window (LUT) center
- window (LUT) width
- zoom factor
- system name
- description of study (if entered)
- progressive number of series in associated study

The information column to the right of the image displays:

- patient's last name and name
- patient code (ID)
- patient's date of birth and sex
- limb district and laterality
- examination date
- study hour
- Institute name
- Department name

The central area displays the orientation of the image with respect to the patient: F refers to feet, H to head or P for posterior and A for anterior. The information regarding the image orientation refer to a patient in a pre-set position according to the DICOM®, standard. This patient is in the supine position, palms facing upwards, as shown in the following figure.


fig. 2.9 - DICOM® , convention for patient positioning

Image layout settings

The image layout can be customized as required, i.e. information can be eliminated or reduced as required.

Proceed as follows:

click with the left mouse button on the corresponding icon on the tool bar alongside: all the information on the image layout disappears, except for the DICOM® position indicators and slice number/total slice number specifications. Furthermore, according to the configured size of the window containing the series, it also displays the name of the sequence, patient's name and last name, region and lateral direction.

Icon functionality is sequential: a second click on the left mouse button on the image returns to the automatic settings of the image layout.

Note that the above is applied exclusively to the current series and not saved in the memory.

Browsing between images

Displaying previous and next images

Inside the current window, all images related to the selected patient can be displayed, by using the previous and next image functions.

The previous image can be displayed as follows:

- 1 METHOD 1, select from the main menu **Image** → **Browse** → **Previous** or click on the relative icon of the tool bar displayed alongside.
- 2 METHOD 2, press ← on the keyboard, or the key ↑ in the case of the format 1 × 1.

The next image can be displayed as follows:

- 1 METHOD 1, select from the main menu **Image** → **Browse** → **Next** or click on the relative icon of the tool bar.
- 2 METHOD 2, press → on the keyboard, or the key ↓ in the case of the format 1 × 1.

In format 1 × 1, the display of the previous or next image involves loading the image that precedes or follows the current image in the exploded series inside the window.

With regard to 2 × 2 and custom formats, the meaning of the meaning of current image needs to be underlined.

☛ Current image refers to the currently active image, highlighted with a yellow border.

In 2 × 2 or custom format, the display of the previous image - starting from any image that is not the image in the upper left corner - means making the previous image (immediately to the left) the current image.

When the image is located in the upper left corner and the previous image needs to be displayed, the latter will be loaded in place of the current image, which will move to the box to the right, and likewise for the other images shown.

In 2×2 or custom format, the display of the next image - starting from any image that is not the image in the lower right corner - means making the next image (immediately to the right) the current image.

When the image is located in the lower right corner and the next image needs to be displayed, the latter will be loaded in place of the current image, which will move to the box to the left, and likewise for the other images shown.

Displaying previous and next series

☛ Images from the previous or next series can be loaded into the active window by following the instructions below.

The previous series can be displayed as follows:

1 METHOD 1, select from the main menu **Image** → **Browse** → **Previous series**.

2 METHOD 2, click on the corresponding icon on the tool bar displayed alongside.

The next series can be displayed as follows:

1 METHOD 1, select from the main menu **Image** → **Browse** → **Next series**.

2 METHOD 2, click on the corresponding icon on the tool bar displayed alongside.

If the images are displayed in 1×1 format, the first image in the previous or next series is loaded.

If the images are displayed in 2×2 or custom format, the first images (according to the format) of the previous or next series are loaded in the current window.

When browsing between series, it is not possible to move to the previous or next study. An acoustic signal warns the user when the first or last series in the study has been reached.

Image properties

A window can be displayed, showing the properties of an image currently open in the viewer.

To do this, if more than one series is open at the same time on the display, first select the required series by clicking on the relative title bar and then select from the main menu:

View → Properties

At this point the **Properties** window shown in the following figure appears, providing the user with the series data.


fig. 2.10 - Properties window

Saving images in non-diagnostic formats

Single images (slices) of each series on the removable media (optical disk or CD/DVD) can be exported onto hard disk as required.

To display images in this way, proceed as follows:

- 1 Insert the removable media (optical disk or CD/DVD) in the relative drive.
- 2 If the workspace to the left of the image viewer is not open, from the main menu select **View → Workspace**.
- 3 In the workspace, click on + next the to icon of the patient displayed. This shows all studies performed on the patient and stored on the optical disk or CD/DVD.
- 4 Click on + next to the icon of the study to be displayed. This shows all the series related to this study and stored on the optical disk or CD/DVD.
- 5 Double click with the left mouse button on the series of images to be exported.

- 6 Move within the series using the arrow keys and/or specific commands on the viewer (see chapter 2 “Viewer”) to move to the image to be exported.
- 7 From the menu, select **File** → **Save as**.
- 8 This opens the “Save as..” window for selection of the folder and name under which to save the image. There is also the option to save in the format jpg or bmp.

Warning

JPEG format always involves compression of data, with loss of information.


fig. 2.11 - Save as

To export several images of the same series, perform this procedure from point 6 to point 8 for each of the images to be exported.

CHAPTER 3

Image processing functions

.....

This chapter discusses the imaging functions available to the user.

All image processing functions described in this chapter are run in the current window, and therefore all images in the current series, regardless of whether the view format is 1×1 , 2×2 or custom.

It is important to note that on selection of one of the functions described below, the mouse pointer remains activated on this function until another one is selected or the same one is deselected by means of the left mouse button (in this case the pointer is no longer activated on any function).

Enlarging images

Sometimes it may be necessary to enlarge a portion of an image to make it more legible; for example, in the case of a small structure of specific interest. If this function is to be applied to the current image, follow the procedures described below:

- 1** METHOD 1, select from the main menu **Image** → **Zoom**, or click on the relative icon of the tool bar displayed alongside.
- 2** METHOD 2, click inside the image with the right mouse button, to display the context menu and select the **Zoom** function.

The zoom function changes the mouse pointer from an arrow into a magnifying glass.

Press and hold the left mouse on the selected image to zoom into the entire image by moving the mouse up, or zoom out by moving the mouse down. The size of the area displayed will remain unchanged in zoom mode.


fig. 3.1 - Image zoom (same image, zoomed)

To reset the image zoom, and return to the original dimensions of all of the images in the current series, proceed as follows:

- 1 METHOD 1, select **Image** → **Reset** → **Zoom/Pan**, from the main menu

- 2 METHOD 2, click on the arrow next to the corresponding icon on the tool bar displayed alongside, to select **Reset Zoom/Pan**.

Changing image brightness and contrast

For an ideal view of images in terms of brightness and contrast, the display table of the image can be modified as required (Look-Up-Table or, abbreviated: LUT). LUT variation modes are available to users and may be activated by means of function keys or the relative menus.

There are three different ways to access the LUT function. The first two are shown below:

- 1 METHOD 1, select from the main menu **Image** → **Lut**, or click on the relative icon of the tool bar displayed alongside.
- 2 METHOD 2, click inside the image with the right mouse button, to display the context menu and select the **Reset Lut** function.

The LUT function changes the mouse pointer from an arrow into a cross. At this point, press the middle mouse button and drag the mouse up or down respectively to increase or decrease brightness, releasing the button when the cursor rests on the value required.

Similarly, by dragging the mouse - with the middle key pressed down - to the right, image contrast can be increased, and by dragging to the left, the contrast is decreased; release the button on display of the required value.

To reset modifications to the LUT and thus restore the images in the current series to their status before adjusting the brightness and contrast, proceed as follows:

1 METHOD 1, select **Image** → **Reset** → **Lut** from the main menu

2 METHOD 2, click on the arrow next to the corresponding icon on the tool bar displayed alongside, to select **Reset Lut**.

Image pan

Image “PAN” refers to moving the image within the frame that holds it, positioned in the view window of the series.

Obviously, if the view format is format 1, “PAN” means moving the image inside the window containing it.

There are two ways to access the PAN function:

1 METHOD 1, select from the main menu **Image** → **Pan**, or click on the relative icon of the tool bar displayed alongside.

2 METHOD 2, click inside the image with the right mouse button, to display the context menu and select the **Pan** function.

The PAN function changes the mouse pointer from an arrow into a hand. By pressing and holding the left mouse button on the selected image, the user can move it by moving the mouse within the window; the image is positioned in the currently occupied area when the button is released. The size of the area displayed will remain unchanged in zoom mode.

To reset the movement of the image carried out and to restore the image of the current series to the position before modification, proceed as follows:

1 METHOD 1, select **Image** → **Reset** → **Zoom/Pan** from the main menu

2 METHOD 2, click on the arrow next to the corresponding icon on the tool bar displayed alongside, to select **Reset Zoom/Pan**.

Changing image orientation

In general, images produced with the standard reconstruction algorithm are already oriented correctly and thus this function is not particularly useful. However, if the orientation of an image needs to be modified, proceed as follows.

☛ To rotate the image on its vertical axis, one of the two following methods can be used:

1 METHOD 1, select from the main menu **Image** → **Rotations/Mirror** → **Mirror**.

- 2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon

(Mirror).

☛ To rotate the image on its horizontal axis, one of two different methods can be used:

- 1 METHOD 1, select from the main menu **Image** → **Rotations/Mirror** → **Flip**

- 2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon

(Flip).

☛ To rotate the image counter-clockwise in the current window at 90 degrees, one of two different methods can be used:

- 1 METHOD 1, select from the main menu **Image** → **Rotations/Mirror** → **Rotate 90° anticlockwise**.

- 2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon

(Rotate 90° anticlockwise).

☛ To rotate the image counter-clockwise in the current window at 180 degrees, one of two different methods can be used:

- 1 METHOD 1, select from the main menu **Image** → **Rotations/Mirror** → **Rotate 180°**.

- 2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon

(Rotate 180°).

☛ There are two different ways to rotate the image in the current window by 90 degree clockwise:

- 1 METHOD 1, select from the main menu **Image** → **Rotations/Mirror** → **Rotate 90° clockwise**.

- 2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon

(Rotate 90° clockwise).


fig. 3.2 - Example of mirror images and rotations

To reset the changes made to image orientation, and restore the images in the current series to their original status, proceed as follows:

1 METHOD 1, select from the main menu **Image** → **Reset** → **Rotations**.

2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon **(Reset Rotations)**.

3 METHOD 3, click on the arrow of the icon displayed alongside to open a context menu and select **Reset Rotations**.

Image Clips

Clips, i.e. cut-outs of a portion of the selected image can be made as required.

The clip is made as follows:

☛ for rectangular clips

1 METHOD 1, select from the main menu **Image** → **Clip** → **Rectangular**.

2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon **(Rectangular)**.

When the left mouse button is pressed and held, the corner of a rectangle is displayed with flexible sides to enable the user to modify the size of the image as required. The size is set on release of the button.

☛ for oval clips

1 METHOD 1, select from the main menu **Image** → **Clip** → **Oval**.

2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon  (**Oval**).

When the left mouse button is pressed and held, the focus of an oval is shown. Its perimeter is flexible and its size varies when the mouse is moved over the image. The size is set on release of the button.

☛ for manual clips

1 METHOD 1, select from the main menu **Image** → **Clip** → **Manual**.

2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon  (**Manual**).

First click the left mouse button at the point to begin drawing the clip, then release the button and use the mouse cursor to outline the first side of the piece to be defined and click again on the left mouse button. Then proceed in the same way with the second side of the piece and so on. The clip is closed automatically on a double click of the mouse, linking the last point to the initial point of the clip with a straight line.

To reset the clipped portion, and restore the images in the current series to their original status, proceed as follows:

1 METHOD 1, select from the main menu **Image** → **Reset** → **Clip**.

2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon  (**Reset Clip**).

3 METHOD 3, click on the arrow of the icon displayed alongside to open a context menu and select **Reset Clip**.

Statistics

A number of functions can be executed on the selected image to calculate or display statistical data regarding a specific area of an image.

In particular, the average and standard deviation of the selected region of interest (ROI) are calculated.

The parameters above are determined by clipping the region of interest according to the procedure below:

☛ for a rectangular ROI

1 METHOD 1, select from the main menu **Image** → **Statistics** → **Rectangular**.

2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon  (**Rectangular**).

When the left mouse button is pressed and held, the corner of a rectangle is displayed with flexible sides to enable the user to modify the size of the image as required. The size is set on release of the button.

The average and standard deviation appear - and vary in real time - next to the rectangle as the size is changed.

☛ with an oval ROI

1 METHOD 1, select from the main menu **Image** → **Statistics** → **Oval**.

2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon  (**Oval**).

When the left mouse button is pressed and held, the focus of an oval is shown. Its perimeter is flexible and its size varies when the mouse is moved over the image. The size is set on release of the button.

The average and standard deviation appear - and vary in real time - next to the oval as the size is changed.

☛ with a manual ROI

1 METHOD 1, select from the main menu **Image** → **Statistics** → **Manual**, or click on the corresponding icon on the tool bar displayed alongside.

2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon  (**Manual**).

Press and hold the left mouse button at the point to begin drawing the ROI, then release the button and move the mouse to outline the first side of the piece to be defined with another click on the left mouse button. Then proceed in the same way with the second side of the piece and so on. The clip is closed automatically on a double click of the mouse, linking the last point to the initial point of the clip with a straight line. The average and standard deviation appear - and vary in real time - from the moment in which you begin drawing the second side of the ROI.

A drawn region of statistical interest can be copied as required. Click the left mouse button on the statistics to be copied - this displays the relative details - and then:

- 1 METHOD 1, select from the main menu **Image** → **Statistics** → **Copy**.
- 2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon  (**Copy**).

At this point a copy of the selected statistics will be displayed with details, and it can be moved to the required position - by dragging it with the mouse, taking care to keep the left mouse button pressed.

To reset the region of interest that has been drawn, and restore the images in the current series to their original status, proceed as follows:

- 1 METHOD 1, select from the main menu **Image** → **Reset** → **Measurements**.
- 2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of **Reset Measurements**.


fig. 3.3 - Example of rectangular statistics

Measurements

Measuring the distance between pairs of points

The distance between multiple pairs of points can be measured as follows:

- 1 METHOD 1, from the main menu, select **Measurements** → **Distances**
- 2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon  (**Distances**).

Move the mouse pointer over the image, where it takes on the shape of a cross, identify the position of the first point and click with the left mouse button on the point. Press and hold the mouse button, and move the mouse pointer to the second point. A flexible segment is outlined during this process. On release of the mouse button, the drawn segment is set and the distance is measured between the two points. This value appears next to the image.

The ends of the segment are highlighted with two small boxes; click on the latter to make further changes as required.

Click on the image again to remove the edit frames.

More than one distance can be displayed between various pairs of points at the same time.

Distances measured this way are approximate. The accuracy increases the closer the pointer is to the magnet isocenter. To obtain measurements with absolute precision, refer to a known example.

fig. 3.4 - Measuring the distance between pairs of points

Measuring the angle between two segments

The angle made by any two segments can be measured as follows:

- 1 METHOD 1, from the main menu, select **Measurements** → **Angle**
- 2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon  (**Angle**).

Move the cursor to the image and it takes on the form of a cross, after which set the starting point of the first segment by clicking on the left mouse button.

Move the mouse to the point for the second end of the segment, which is also as the vertex of the angle to measure; click to set this point.

Move the mouse to set, with a click, the second end of the second segment.

The ends of the segments and the vertex of the angle are highlighted with frames; if required, click on these frames to make further changes.

Click on the image again to remove the edit frames.

More than one angle of various segment pairs can be displayed at the same time.

The value of the acute angle between two segments appears - and varies in real time - from the moment the user starts to draw the second segment.

The angle value measured in this way is approximate. The accuracy increases the closer the pointer is to the magnet isocenter. To obtain measurements with absolute precision, refer to a known example.

Inserting text and direction arrows

The following can be inserted in images:

☛ text

- 1 METHOD 1, from the main menu, select **Measurements** → **Text**
- 2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon  (**Text**).

Then click with the left mouse button on the point in the image for text entry and enter the text via the keyboard.

The text is framed in a rectangle with ends highlighted with frames; click on the latter to edit the text as required.

Click on the image again to remove the edit frames.

☛ arrows

- 1 METHOD 1, from the main menu, select **Measurements** → **Arrow**
- 2 METHOD 2 click on the arrow of the icon displayed alongside to open the icon context menu, for selection of the icon  (**Arrow**).

Then click with the left mouse button on the point for arrow insertion; press and hold the button and drag the arrow to the required position.

The arrow is contained in a rectangle with ends highlighted with frames; click on the latter to make further modifications to the arrow as required.

Click on the image again to remove the edit frames.

fig. 3.5 - Inserting text and arrows in an image

Deleting measurements taken

To reset measurements of an image and restore the images in the current series to their original status, proceed as follows:

- 1 METHOD 1, select the measurement to be deleted, followed by **Measurement** → **Cancel** from the main menu

- 2 METHOD 2, select the measurement to be deleted, then click on the arrow of the icon displayed alongside to open the icon context menu, and select the icon  (**Cancel**).

- 3 METHOD 3, select the measurement to be deleted and press **CANC** on the keyboard.
- 4 METHOD 4 select **Image** → **Reset** → **Measurements** from the main menu to cancel all measurements.

- 5 METHOD 5, click on the arrow of the icon displayed alongside to open the icon context menu, and select **Reset Measurement** to delete all measurements made

- 6 METHOD 6, click on the arrow of the icon displayed alongside to open the icon context menu, and select **Reset Measurement** to delete all measurements made

Deleting all changes made to an image

If various image processing functions have been carried out on an image (LUT change, insertion of measurements, statistics, rotations, etc) these can all be deleted, if required, to return to the original image.

- 1 METHOD 1, click on the arrow of the icon shown alongside, to open the sub-menu and select **Reset All**

- 2 METHOD 2, click on the icon displayed alongside (on the icon and not the arrow)
- 3 METHOD 4, select from the main menu **Image** → **Reset** → **All**.

CHAPTER 4

User Interface management

.....

Showing / hiding tool bars

The user interface can be configured according to user requirements or habits.

The **Viewer** toolbar is always active and cannot be hidden, with the following icons:

There is also the option to show / hide the status bar as described below.
From the main menu, select **View** → **Status Bar**
once to show it if hidden and vice versa.

Window management

When more than one window is opened at a time, the system, by default, places the relative windows side by side inside the image viewer: the second window opened is placed horizontally next to the first, which is redimensioned, after which if other series are opened, the display follows the tile display as adopted in the Windows®, operating system, and shown in the following figure.

fig. 4.1 - "Tile" arrangement of windows

If the user changes the arrangement of a window or its size, the viewer no longer automatically manages window positioning.

The windows can be re-arranged by default at any time:

select **Window** → **Auto Arrange**, from the main menu or click on the relevant icon on the toolbar.

The windows can also be positioned in cascade:

select **Window** → **Cascade**, from the main menu or click on the relevant icon on the toolbar.

fig. 4.2 - Displaying different series in cascade

To return to the automatic tile arrangement of windows:

select **Window** → **Auto Arrange**, from the main menu or click on the relevant icon on the toolbar.

With both the tile and cascade display, any window can be maximized to full screen as required.

To maximize a window, click with the left mouse button on the window's maximize/minimize button, or double click on the title bar.

To minimize a window, click with the left mouse button on the window's maximize/minimize button.

Closing a window

A window can be closed in two ways. First ensure that the window to be displayed is the current one on screen.

If the window to be closed is not the current one on display, click with the left mouse button on the blue bar placed at the top of the required window, and proceed as follows:

- 1 METHOD 1, from the main menu, select **File** → **Close**
- 2 METHOD 2, click with the left mouse button on the symbol "X" at the top right corner of the window.

About VLite

From the main menu, select ? → **About VLite** from the main menu, a window is displayed with the following information:

- the software revision, with the relevant date of release


fig. 4.3 - **About VLite** window

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